

PLASTER OF PARIS BANDAGES PLATROFIX



INDICATIONS :

- Fractures
- Sprains
- Dislocations
- Orthopedic Surgeries...

CHARACTERISTICS :

- Quick and regular moistening
- Water temperature :20-25°C
- Creamy plaster
- Frequent smoothing
- Good modeling properties
- Setting time : 3 to 6 min
- Short drying time
- Colour: White

DIMENSIONS :

- 8 cm * 3 m
- 10 cm * 3 m
- 15 cm * 3 m
- 20 cm * 3 m

OPTIONS (On request) :

- Choice of colors
- Choice of Fragrances



- ### Advantages of the PLATROFIX Bandage
- Platrofix is made of 100% cotton gauze, impregnated with plaster and binder, rolled on a plastic tube, non-sterile and for single use only.
 - Platrofix is designed for quick fixation of plaster with reduced immersion and excellent casting ability with low loss of plaster. Once completely hardened, becomes very resistant although consumes little bandages.
 - Platrofix is widely used due its effectiveness. It is effective in immobilizing fractures and promoting bone healing. As for its effectiveness, it is no longer to be proved.

